

WHAT IS ACCELEDENT AURA?

AcceleDent® Aura is a simple-to-use, hands-free device designed for faster orthodontic treatment—with only 20 minutes daily use.



» It's Fast.

Clinically proven to move teeth up to 50% faster in a US clinical trial.¹

» It's Safe.

Cleared by the United States Food and Drug Administration and available only by prescription from your orthodontist.¹

» It's Gentle.

AcceleDent patients have reported reduced discomfort with their orthodontics.¹

Frequently Asked Questions

How does AcceleDent Aura work?

By inserting the mouthpiece fitted around your existing orthodontics and wearing the activator for 20 minutes everyday, AcceleDent Aura can speed up tooth movement through the use of SoftPulse Technology™. Faster tooth movement may decrease the duration of your orthodontic treatment and may help to make your orthodontics more comfortable.¹

Is AcceleDent Aura safe?

Yes, AcceleDent has been proven safe in clinical trials and is cleared by the FDA as a Class II medical device.¹

What is SoftPulse Technology™ and why is it important?

Part of the way that orthodontics works is by changing or remodeling the bones surrounding your teeth. As this remodeling process is accelerated, your teeth move more quickly. That's where AcceleDent Aura comes in. By using tiny vibrations or micropulses, SoftPulse Technology™ speeds up bone remodeling and tooth movement.

TALK TO YOUR ORTHODONTIST ABOUT ACCELEDENT:



Visit www.aceledent.com for more information
or call 1-866-866-4919.



AcceleDent is an orthodontic accessory intended for use during orthodontic treatment. Refer to the Directions for Use with AcceleDent Aura or aceledent.com for warnings and precautions prior to using this product. Federal (U.S.) law restricts this device to sale by or on the order of a licensed doctor. Rx use only for a single patient.

Made in the USA

©2013 OrthoAccel® Technologies, Inc.

AcceleDent® is a registered trademark of OrthoAccel® Technologies, Inc.

Oral-B® is a registered trademark of P&G.



Reference: 1. Data on File.

L61000 Rev.2, JUNE 2013

FASTER ORTHODONTIC TREATMENT



NEW
AcceleDent®
aura




» SoftPulse Technology micropulses are gentle and exert up to 8x less force than a power toothbrush.




■ AcceleDent Aura

vs.



■ Oral-B® Sonic Complete power toothbrush



» Your AcceleDent Aura comes with four ergonomically friendly components:

ACTIVATOR

The heart of AcceleDent Aura, the small, lightweight Activator generates gentle micro-pulses. It includes a USB interface which can be plugged directly into a computer to view your usage history via the FastTrac Report.



MOUTHPIECE

Chosen specifically for you by your orthodontist, the Mouthpiece provides a comfortable fit and snaps easily on and off the Activator for transport and cleaning.



USB EXTENSION CABLE with POWER ADAPTER

Allows for easy, convenient charging.



TRAVEL CASE

The Travel Case keeps the components safe, clean and secure.



ACTUAL PATIENT RESULT

Patient: Fred C.

AcceleDent Treatment Time: 9.5 months

Estimated Treatment Time without AcceleDent: 18 months

Robert T. Rudman, DDS, MS

Cherry Creek Smiles, Denver, CO